en TO REPLACE SIEVE BEDS

⚠ CAUTION– Always work in an ESD (electrostatic discharge) safe environment when assembling or disassembling this device.

1. Ensure that all power is removed from the unit including battery.
2. Lay the unit on its side with the USB port facing upward.
3. Using a T-10 torx screwdriver, remove all 6 screws from the side case and remove the case cover from unit. (Fig. 1)
4. Carefully lift the front display cover plate up and out of case. Place it beside unit. (Fig. 2)
5. Using a T-10 torx screwdriver, loosen the 2 screws at the top of the valve manifold of the sieve bed module. (Fig. 3)
6. Lift sieve bed module up several inches out of the retaining catches and remove the valve manifold from the top of the sieve beds. (Fig. 4A & 4B)
7. Remove and discard the original sieve beds and all four rubber O-Rings located at the top and bottom of each sieve bed. (Fig. 5)
 - ⓘ **NOTE**– The O-Rings may remain on the valve manifold or in the recesses of the accumulator tank.
8. Uncap two new sieve beds and seat the bottom O-Rings (larger O-Rings) on the bottom of each sieve bed, and the top O-Rings on the aluminum manifold (smaller O-Rings). (Fig. 6)
 - ⓘ **NOTE**– 4 O-Rings are included with each sieve bed package.
9. Place sieve beds over accumulator tank assembly ports. Then rotate beds until the bed is fully seated into position and so the serial number is viewable. Be careful to prevent pinching the O-Rings during assembly. (Fig. 7A & 7B)
10. Gently lift the sieve beds and accumulator tank upwards to re-insert the manifold onto the top of the sieve beds. Rotate the sieve beds slightly while inserting to prevent pinching of the O-Rings.
11. Re-install the 2 manifold screws (torque to 6-8 in-lbs). (Fig. 8)
12. Lift the front display cover plate back into the slots of the side case. (Fig. 9)
13. Ensure that no hoses or wiring harnesses are pinched and place the side cover back onto the unit.
14. Re-install the 6 case screws (torque to 7 in-lbs).
15. Re-install the battery into the unit. (Fig. 10)
 - ⓘ **NOTE**– iGo2 automatically increases compressor speed when sieve beds are at the end of their life cycle. The increased compressor speed will continue after the sieve beds are replaced and is normal. This does not cause any safety or performance risk. The compressor speed will normalize after a short period of operation of the unit. After installing new sieve beds, the iGo2 should be operated without the cannula connected for 20 minutes at each setting.
16. If the unit still does not function properly after reassembly with new sieve beds:
 - a. Remove the new sieve beds and reinstall cap and plug. To preserve sieve beds for an extended period, place them in a sealed bag after capping. Note: Exposure to ambient air can damage sieve and negatively affect POC performance.
 - b. Follow return procedure to send the unit to Drive DeVilbiss Healthcare or an Authorized Service Center for servicing; the defective sieve beds should be returned with the unit.

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