

# FreeStyle® Comfort®

The perfect fit



## Key features and benefits



### Optimal Comfort

Let your patients experience the comfort of the body-skimming, curved silhouette design. The FreeStyle Comfort portable oxygen concentrator is lightweight and compact.



### Smart O2 delivery

Unique clinical features including UltraSense®, autoDOSE®, and myCAIRE™ help to ensure your patients will get the oxygen they need, when they need it.



### Proper O2 Saturation

With 5 settings and up to 1050 mL of oxygen per minute, the FreeStyle Comfort serves a wide range of patient needs. The unit delivers up to 95% pure oxygen to further ensure proper saturation.



### Proven Quality

Robust design features, including the sieve bed technology and critical components, allow for minimum maintenance and long life. CAIRE offers one of the best warranties available, with a 3-year warranty on the device and sieve beds.

## Specifications



### FLOW RATE

Pulse Dose 1-5  
O2 Output Settings:  
Setting 1-210 mL/min  
Setting 2-420 mL/min  
Setting 3-630 mL/min  
Setting 4-840 mL/min  
Setting 5-1050 mL/min

### DIMENSIONS

10.0 x 7.3 x 3.1 in  
(25.4 x 18.5 x 7.9 cm)  
8-cell battery pack  
11.0 x 7.3 x 3.1 in  
(27.9 x 18.5 x 7.9 cm)  
16-cell battery pack

### WEIGHT

5 lbs (2.3 kg)  
8-cell battery pack  
6 lbs (2.7 kg)  
16-cell battery pack

### SOUND LEVEL

39.9 dB(A) @ Setting 2

### BATTERY DURATION

8-cell battery  
up to 8 hours  
16-cell battery  
up to 16 hours

### POWER

AC Power: 100-240  
VAC, 50-60Hz  
DC Power: 11-18 VDC  
(10 max amp)  
Battery: Lithium Ion

### WARRANTY

3 years  
(Includes sieve beds)  
1 year on battery and  
power supplies

# FreeStyle® Comfort®

Open the door to a bigger,  
more enjoyable life

**CARRYING BAG  
(Holds FreeStyle Comfort only)**

**Item Number**  
MI396-1



**AC INPUT CABLE  
(Regional Options Available)**

**Item Numbers**  
4997-SEQ (Type A Plug)\*  
4998-SEQ (Type C Plug)\*



**CARRY-ALL ACCESSORY BAG  
(Holds Unit and Accessories)**

**Item Number**  
MI320-2

**Replacement Strap  
Item Number**  
MI397-1



**BATTERY PACK**

**8-Cell Item Number**  
BT034-2S



**16-Cell Item Number**  
BT035-1S



**DESKTOP BATTERY CHARGER  
(Battery Not Included)**

**Item Numbers**  
BT036-1  
(Type B Plug)\*



BT036-2  
(Type E Plug)\*\*



**BACKPACK**

**Item Number**  
MI459-1



**VEHICLE POWER CORD  
(12 VDC Outlet)**

**Item Number**  
CD041-1



**TRAVEL ADAPTER**

**Item Number** 5591-SEQ



**POWER SUPPLY WITH  
AC INPUT CABLE**

**Item Numbers**  
PW036-1S  
(Type A Plug)\*



PW036-2S  
(Type C Plug)\*

**POWER SUPPLY (120W 15V)**

**Item Number**  
PW036-1



**GROSS PARTICLE FILTERS**

<b>Right Filter Item Number</b> FI226-1	<b>Left Filter Item Number</b> FI227-1
	



\* See <https://www.iec.ch/worldplugs/map.htm#> for details on plug types.

\*\* Purchase 5591-SEQ along with this item to properly adapt the plug for IEC type I power sockets.

Find out more about FreeStyle Comfort and accessory products at **CAIREinc.com**

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2025 CAIRE Inc. All Rights Reserved. ML-CONC0159 E