



Oxlife **INDEPENDENCE**<sup>®</sup>

# Proven best-in-class reliability, year after year.



Meets the majority of long-term oxygen therapy patient requirements



Continuous flow up to 3 LPM with large pulse flow



Long-life, dual battery operation of up to six hours



Proven 24/7 use, year after year - backed by a five-year warranty



Use as both a POC and stationary unit




No tanks to refill or the limitations of ongoing tank deliveries

## Designed for non-delivery

- The foundation device for your oxygen fleet
- Where intelligent economics meets enhanced patient care
- Quick set-up for hospital discharge



- Device initiated performance monitoring
- Clinical monitoring and compliance support
- Inventory management with Geofence



AMERICAN INNOVATION:  
OWNERSHIP  
ENGINEERING  
MANUFACTURING  
CUSTOMER SUPPORT

ENERGY SMART TECHNOLOGY<sup>®</sup>

# Oxlife INDEPENDENCE<sup>®</sup>

## PRODUCT SPECIFICATIONS<sup>1</sup>

<b>Dimensions</b>	Height: 20.2" • Width: 10.9" • Depth: 9.45"
<b>Weight<sup>2</sup></b>	Device: 16.7 lbs
<b>Oxygen Flow Rate</b> (Pulse 15-40 BPM)	Continuous: 0.5 - 3 LPM in 0.5 increments Pulse: 0.5 - 6.0 in 0.5 increments (8 - 96 ML)
<b>Battery Duration<sup>3</sup></b>	5.75 hours @ 2 pulse at 20 bp 2.5 hours @ 2 LPM continuous
<b>Battery Recharge Time<sup>4</sup></b>	2.5 hours
<b>Oxygen Concentration</b>	87% - 95% at sea level at all settings
<b>AC Power</b>	Universal Input: 100 - 240V, 50/60Hz 50 - 60 Hz • 2.0A maximum
<b>DC Power</b>	12 - 15 VDC • 15 amps
<b>Operating Temperature</b>	50°F (10°C) to 104°F (40°C)
<b>Altitude Operating Range</b>	0 - 13,123 feet (0 - 4,000 meters)
<b>Transportation/Storage Humidity</b>	0 - 95% non-condensing
<b>Connectivity - DNA Technology</b>	O <sub>2</sub> Concepts <sup>®</sup> proprietary, device-initiated, asset and clinical management platform utilizing the Verizon <sup>®</sup> cellular network

<sup>1</sup> O<sub>2</sub> Concepts<sup>®</sup> reserves the right to revise specifications at anytime without notice

<sup>2</sup> Production weight excludes cart kit

<sup>3</sup> Battery duration presented for two-battery system

<sup>4</sup> Battery charge time increases if device is operating while charging



[o2-concepts.com](https://o2-concepts.com) | 877.867.4008

© 2023 O<sub>2</sub> Concepts, LLC; all rights reserved; O<sub>2</sub> Concepts<sup>®</sup> and Oxlife INDEPENDENCE<sup>®</sup> are registered trademarks of O<sub>2</sub> Concepts, LLC and may not be reproduced without express consent.

