

# P2-E7 PORTABLE OXYGEN CONCENTRATOR



ITEM #	DESCRIPTION
P2-E7	7 Setting Portable Oxygen Concentrator

- Compact & lightweight.
- High efficiency, low noise. Excellent heat dissipation.
- 93% O<sub>2</sub> concentration.
- High adsorption sieve bed.
- Includes AC/DC Car Charger.
- Conforms to all applicable FAA compliant criteria for POC carriage and use on board aircraft.

## FEATURES

Size	8.7" x 3.4" x 6.3"
Weight	4.37lbs
Warm up time	2 minutes
Settings	Seven Flow Settings (1-7)
Maximum Oxygen Output	1400 ml/min
Oxygen Concentrator	87%-96%
AC Power Supply:	100-240v, 50-60hz (auto sensing to allow for worldwide use)
DC Power Supply:	Input: 11-16VDC, Output: 19V 6.3a
Rechargeable Battery	14.4VDC
Battery Rated Capacity:	6.7Ah
Display Screen	2.8"
Battery	6700mAh, duration up to 6 hours on setting 1. Approx. 4 hours charge time with AC.
Operation	Simple control functions and easy-to-read LCD display

Setting	Battery Life
1	6 h
2	5.2 h
3	3.6 h
4	2.8 h
5	2.3 h
6	1.9 h
7	1.4 h

## WARRANTY

Three years, excluding battery and sieve bed

V10: 9/16/24