

P2-S4 PORTABLE OXYGEN CONCENTRATOR



ITEM #	DESCRIPTION
P2-S4	Portable Oxygen Concentrator

- Compact & Lightweight.
- High Efficiency, Low Noise. Excellent Heat Dissipation.
- High Adsorption Sieve Bed.
- Long Lasting Lithium-Ion Battery - Up to 7.2 Hours Battery Life.
- Single Battery Standard. Available Optional Double and Triple Battery.
- Includes AC Adapter & DC Adapter (Car Charger).
- Conforms to all applicable FAA acceptance criteria for POC carriage and use on board aircraft.

SPECIFICATIONS:

Size	6.2" x 3.4" D x 6.7" H (with Single Battery) 6.2"L x 3.4"W x 7.5"H (with Double Battery) 6.2"L x 3.4"W x 8.3"H (with Triple Battery)
Weight	3.3 lbs (with Single Battery) 4 lbs (with Double Battery) 4.4 lbs (with Triple Battery)
Warm-up Time	2 Minutes
Settings	Four Flow Settings (1-4)
Oxygen Concentration	87%-96%
AC Power Supply	100-240v, 50-60hz
DC Power Supply	1-16V DC Output: 19V 6.3A Voltage: 14.4V
Rechargeable Battery	14.4V DC 3.35A MAX Single Battery Rated capacity: 3,400 mAh Double Battery Rated capacity: 6,800 mAh Triple Battery Rated capacity: 10,200 mAh Duration up to 2.5 hours on setting 1 Approx. 2 hours charge time with AC
Display Screen	2"
Operation	Simple Control Functions and Easy-to-Read LCD Display
Use	Designed for 24/7 use

BATTERY DURATION

Single Battery

- Setting 1: 2.5h
- Setting 2: 1.7h
- Setting 3: 1.2h
- Setting 4: 1h

Double Battery

- Setting 1: 4.8
- Setting 2: 3.5
- Setting 3: 2.5h
- Setting 4: 2h

Triple Battery

- Setting 1: 7.2h
- Setting 2: 5h
- Setting 3: 3.5h
- Setting 4: 3h

WARRANTY

Three Years. Sieve Bed: Two Years. Battery: One Year

V6: 10/3/24